MAYO CLINIC LABORATORIES

Test Definition: FBACS

Baclofen, Serum

Overview

Method Name
High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Serum Red

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Collection Container/Tube: Red top (Serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collections Instructions:
1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
Serum: 0.4 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen	120 days	

Clinical & Interpretive

Clinical Information

Refer to www.nmslabs.com/test-catalog

Reference Values
Reporting limit determined each analysis.
None Detected

Interpretation
Serum concentrations required for therapeutic effects reportedly range from 0.08-0.40 mcg/mL.

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
5 to 11 days

Performing Laboratory Location
NMS Labs

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information
80369

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBACS	Baclofen, Serum	9353-4

Result ID	Test Result Name	Result LOINC® Value

Z3532	Baclofen	9353-4