MAYO CLINIC LABORATORIES

Test Definition: FHCC1

Hazelnut Component rCor a 1

Overview

Method Name
ImmunoCAP FEIA

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red or SST
Specimen Volume: 0.5 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume
0.5 mL

Reject Due To

Thawing:	Warm OK; Cold OK
Gross lipemia:	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values
<0.10 kU/L

Interpretation

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10 - 0.34	Equivocal/Borderline
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	High Positive
4	17.50 - 49.99	Very High Positive
5	50.00 - 99.99	Very High Positive
6	>99.99	Very High Positive

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
3 to 6 days

Performing Laboratory Location
Eurofins Viracor

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test has been cleared or approved for diagnostic use by the U.S. Food and Drug Administration.

CPT Code Information
86008

LOINC® Information

MAYO CLINIC LABORATORIES

Test Definition: FHCC1

Hazelnut Component rCor a 1

Test ID	Test Order Name	Order LOINC® Value
FHCC1	Hazelnut Component rCor a 1	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5658	Hazelnut Component rCor a 1	69421-6
Z5659	Class	103076-6