MAYO CLINIC LABORATORIES

Test Definition: FFRFT

Anti-Synthetase Profile

Overview

Method Name
RIPA Gel Radiography; Enzyme-Linked Immunosorbent Assay (ELISA)

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube: Red top or serum gel
Submission Container/Tube: Plastic vial
Specimen Volume: 3 mL Serum
Collection Instructions:
1. Within 1 hour of collection, centrifuge and aliquot 3 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
Serum: 1 mL; Note: This volume does not allow for repeat testing.

Reject Due To

Gross hemolysis:	Reject
Lipemia:	Reject
Icterus	Reject
Anything other than serum	Reject
Bacterial contamination	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	7 days	
	Frozen	60 days	

Clinical & Interpretive

Clinical Information
Refer to www.labcorp.com/test-menu

Reference Values
Anti-PL-7 Ab, Anti-PL-12 Ab, Anti-EJ Ab, Anti-OJ Ab: Negative

Anti-Jo-1 Ab:
Negative: <20 units
Weak Positive: 20-39 units
Moderate Positive: 40-80 units
Strong Positive: >80 units

Performance

PDF Report
No

Day(s) Performed
Monday through Friday

Report Available
16 to 25 days

Performing Laboratory Location
Esoterix Endocrinology

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification
This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

CPT Code Information
83516 x 4

86235

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFRFT	Anti-Synthetase Profile	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z5782	Anti-Jo-1 Ab	35333-4
Z5783	Anti-PL-7 Ab	33772-5
Z5837	Anti-PL-12 Ab	33771-7
Z5823	Anti-EJ Ab	45149-2
Z5824	Anti-OJ Ab	45152-6