MAYO CLINIC LABORATORIES

Test Definition: FAS1A

Anti-SAE1 Ab, IgG

Overview

Method Name
Enzyme-Linked Immunosorbent Assay (ELISA)

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Collection Container/Tube: Serum gel or red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Within 1 hour of collection, centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume
Serum: 0.5 mL Note: This volume does NOT allow for repeat testing.

Reject Due To

Gross hemolysis:	Reject
Thawing:	Cold OK; Warm OK
Gross lipemia:	Reject
Gross icterus	Reject
Bacterial contamination	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Ambient	7 days	
	Frozen	60 days	

Test Definition: FAS1A

Anti-SAE1 Ab, IgG

Clinical & Interpretive

Clinical Information

Refer to www.labcorp.com/test-menu/

Reference Values

<20
Negative: <20 units
Weak Positive: 20-39 units
Moderate Positive: 40-80 units
Strong Positive: >80 units

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

11 to 14 days

Performing Laboratory Location

Esoterix Endocrinology

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

CPT Code Information

83516

LOINC® Information

MAYO CLINIC LABORATORIES

Test Definition: FAS1A

Anti-SAE1 Ab, IgG

Test ID	Test Order Name	Order LOINC® Value
FAS1A	Anti-SAE1 Ab, IgG	101220-2

Result ID	Test Result Name	Result LOINC® Value
FAS1A	Anti-SAE1 Ab, IgG	101220-2