MAYO CLINIC LABORATORIES

Test Definition: FLISD

Lisdexamfetamine as Metabolite, Urine

Overview

Method Name
High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Urine

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Container/Tube: Plastic preservative-free urine container
Specimen Volume: 1 mL
Collection Instructions:
1. Send 1 mL random urine specimen in a plastic preservative-free urine container.
2. Ship refrigerate.

Specimen Minimum Volume
0.2 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	180 days	
	Ambient	7 days	
	Frozen	180 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog.

Reference Values

Reporting limit determined each analysis.

None detected

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

9 to 13 days

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80324

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FLISD	Lisdexamfetamine as Metabolite, U	19346-6

Result ID	Test Result Name	Result LOINC® Value
Z5825	Amphetamine	19346-6