MAYO CLINIC LABORATORIES

# Test Definition: FFIG2

IGF-2

## Overview

### Method Name
ELISA

### NY State Available
Yes

## Specimen

### Specimen Type
Serum

### Specimen Required
**Collection Container/Tube:** Serum gel or red top
**Submission Container/tube**: Plastic vial
**Specimen Volume**: 0.5 mL Serum
**Collection Instructions**:
1. Within 1 hour of collection, centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Freeze immediately after separation and send frozen.

### Specimen Minimum Volume
Serum: 0.3 mL; Note: Minimum volume does not allow for repeat analysis.

### Reject Due To

| | |
|---|---|
| Specimen not serum | Reject |

### Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---|---|---|---|
| Serum | Frozen (preferred) | 14 days | |
| | Ambient | 14 days | |
| | Refrigerated | 14 days | |

## Clinical & Interpretive

### Clinical Information
Refer to www.labcorp.com/test-menu

### Reference Values
Prepubertal: 258-882 ng/mL
Pubertal: 273-892 ng/mL
Adults: 333-967 ng/mL

## Performance

### PDF Report
No

### Day(s) Performed
Varies

### Report Available
11 to 17 days

### Performing Laboratory Location
Esoterix Endocrinology

## Fees & Codes

### Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

### Test Classification
This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

### CPT Code Information
83520

### LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---|---|---|
| FFIG2 | IGF-2 | Not Provided |

| Result ID | Test Result Name | Result LOINC® Value |
|---|---|---|
| FFIG2 | IGF-2 | Not Provided |