MAYO CLINIC LABORATORIES

Test Definition: FNAD

Nadolol, Serum/Plasma

Overview

Method Name
High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
*****Must submit one specimen per order. Specimens cannot be shared between multiple orders.*****

Submit only 1 of the following specimens:

Specimen Type: Serum
Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Specimen Volume: 1 mL Serum
Collection Instructions:
1. Promptly centrifuge and aliquot serum into a plastic vial.
2. Label specimen as serum.
3. Send refrigerated.

Specimen Type: Plasma
Collection Container/Tube: Lavender top or pink top (EDTA) (plasma gel/PST are **not acceptable**)
Specimen volume: 1 mL Plasma
Collection Instructions:
1. Promptly centrifuge and aliquot plasma into a plastic vial.
2. Label specimen as plasma.
3. Send refrigerated.

Specimen Minimum Volume
0.4 mL

Reject Due To
All specimens will be evaluated by the processing and performing laboratories for test suitability.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	30 days	
	Ambient	30 days	
	Frozen	30 days	

Clinical & Interpretive

Clinical Information
Refer to www.nmslabs.com/test-catalog.

Reference Values
Reporting limit determined each analysis
None detected

Mean steady-state plasma levels following a daily regimen:
80 mg: 26-36 ng/mL
160 mg: 52-74 ng/mL
320 mg: 154-191 ng mL

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
8 to 10 days

Performing Laboratory Location
NMS Labs

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

Test Definition: FNAD

Nadolol, Serum/Plasma

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80299

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FNAD	Nadolol	12406-5

Result ID	Test Result Name	Result LOINC® Value
Z1430	Nadolol	12406-5
Z1867	Reporting Limit	19147-8