MAYO CLINIC LABORATORIES

Test Definition: FHSAG

MVista Histoplasma Ag Quantitative, Spinal Fluid

Overview

Method Name
Quantitative Sandwich Enzyme Immunoassay (EIA)

NY State Available
Yes

Specimen

Specimen Type
CSF

Specimen Required
Specimen Type: Spinal Fluid
Sources: CSF
Container/Tube: Sterile container
Specimen Volume: 0.8 mL
Collection Instructions: Collect 0.8 mL of spinal fluid (CSF). Ship refrigerated, 0.8 mL of spinal fluid. Send specimen in a plastic, screw-capped vial refrigerated.

Specimen Minimum Volume
0.8 mL

Reject Due To

Other	Specimen that is too viscous to pipette. Tissue, biopsy, sputum, bronchial brush, tracheal aspirate, FNA, bone marrow aspirate, stool or samples in transport media, fixative or Isolator tubes

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Refrigerated (preferred)	14 days	
	Ambient	14 days	
	Frozen		

Clinical & Interpretive

Reference Values

Reference interval: None Detected

Reportable Range: Positive Results reported in ng/mL from 0.20 ng/mL to 20.00 ng/mL.

Positive Results above 20.00 ng/mL are reported as “Above the Limit of Quantification”.

Cautions

Cross-reactions are seen with blastomycosis, paracoccidioidomycosis, penicilliosis, less frequently in coccidioidomycosis, rarely in aspergillosis and possibly sporotrichosis.
Sputolysin, sodium hydroxide and potassium hydroxide treatment degrade the analyte detected in the assay.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 5 days

Performing Laboratory Location

MiraVista Diagnostics

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

Test Classification

This test was developed and its performance characteristics determined by MiraVista Diagnostics. It has not been cleared or approved by the FDA; however, FDA clearance or approval is not currently required for clinical use. The results are not intended to be used as the sole means for clinical diagnosis or patient management decisions.

Test Definition: FHSAG

MVista Histoplasma Ag Quantitative, Spinal Fluid

CPT Code Information

87385

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FHSAG	MVista Histoplasma Ag, CSF	51754-0

Result ID	Test Result Name	Result LOINC® Value
Z1722	Result:	51754-0
Z1034	Interpretation	Not Provided