MAYO CLINIC LABORATORIES

Test Definition: FBERY

Beryllium, Blood

Overview

Method Name
Inductively Coupled Plasma/Mass Spectrometry (ICP-MS)

NY State Available
Yes

Specimen

Specimen Type
Whole Blood EDTA - Metal Free (ERB)

Specimen Required
Draw blood in a metal free, royal blue-top with EDTA tube(s). Send 2 mL of EDTA whole blood refrigerated.

Specimen Minimum Volume
0.3 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA - Metal Free (ERB)	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	180 days	

Clinical & Interpretive

Reference Values
Reference Range: <1.0 ng/mL

Test Definition: FBERY

Beryllium, Blood

Performance

PDF Report
No

Day(s) Performed
Monday through Sunday

Report Available
3 to 7 days

Specimen Retention Time
2 weeks

Performing Laboratory Location
Medtox Laboratories, Inc.

Fees & Codes

Fees
- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information
83018

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FBERY	Beryllium	8159-6

Result ID	Test Result Name	Result LOINC® Value
Z1177	Beryllium	8159-6