MAYO CLINIC LABORATORIES

Test Definition: PNZN

Perphenazine, Serum or Plasma

Overview

Method Name
Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
Submit only 1 of the following specimens:

Specimen Type: Plasma
Supplies: Amber Frosted Tube, 5 mL (T915)
Collection Container/Tube: Green-top (sodium heparin) (Plasma gel/PST are **not** acceptable)
Submission Container/Tube: Amber plastic vial
Specimen Volume: 3 mL
Collection Instructions:
1. Draw blood in a green-top (sodium heparin) tube.
2. Centrifuge and aliquot 3 mL plasma into an amber vial to protect from light and refrigerate immediately.
3. Ship refrigerated.

Specimen Type: Serum
Supplies: Amber Frosted Tube, 5 mL (T915)
Collection Container/Tube: Red top (Serum gel/SST are **no**t acceptable)
Submission Container/Tube: Amber Plastic vial
Specimen Volume: 3 mL
Collection Instructions:
1. Draw blood in a plain red-top tube.
2. Centrifuge and aliquot 3 mL of serum into an amber vial to protect from light and refrigerate immediately.
3. Ship refrigerated.

Specimen Minimum Volume
0.5 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	LIGHT PROTECTED
	Ambient	72 hours	LIGHT PROTECTED
	Frozen	180 days	LIGHT PROTECTED

Clinical & Interpretive

Reference Values

Reference Range: 5.0 - 30.0 ng/mL

Low-dose therapeutic range for Perphenazine: 0.5 - 2.5 ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

7 to 9 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to Test Prices for detailed fee information.
- Clients without access to Test Prices can contact Customer Service 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

CPT Code Information

80299

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
PNZN	Perphenazine (Trilafon)	3927-1

Result ID	Test Result Name	Result LOINC® Value
Z1052	Perphenazine	3927-1