


## CYP2D6 Comprehensive Cascade Testing Algorithm\*





<sup>\*</sup>Additional testing is performed to derive as accurate a phenotype as possible.

- 1. No call on an assay for a single nucleotide variant
- 2. Genotype pattern does not fit known star allele
- 3. Duplication involving alleles with different levels of activity
- 4. Genotype/phenotype discord
- 5. Laboratory director request (CYP2D6 real time PCR and CNV assays must be completed first)
- 6. Ambiguous genotype leading to ambiguous phenotype.

<sup>&</sup>lt;sup>A</sup>Always performed

<sup>&</sup>lt;sup>B</sup>Additional testing may be performed (only approximately 3% of individuals may need additional testing) for the following reasons:

<sup>&</sup>lt;sup>c</sup> Tier 2 testing performed on blood specimens only: The number of sequencing tests needed to determine the phenotype will vary depending on the tier 1 result. Due to lower concentration of DNA yielded from saliva, testing cannot proceed to tier 2 sequencing and will stop after tier 1 testing is complete.