



BAP1Z

PATIENT NAME TESTRNV, IMPLEMENTATION				ORDER NUMBER L002000200
PATIENT ID X100366871	DATE OF BIRTH 02/19/1968	AGE 53 Y	SEX Female	REQUESTED BY CLIENT TESTING
COLLECTED 9/2/2021, 7:00 AM	RECEIVED 9/3/2021, 8:09 AM	REPORTED 9/3/2021, 11:56 AM		
The collected, received, and reported dates and times on the report are in the time zone of the performing location. 7028846 MCL RochesterCampus Rochester MN 55901				CLIENT ORDER NUMBER X100366871 CLIENT MRN 321

TEST DESCRIPTION
Evaluation of the BAP1 gene associated with BAP1-tumor predisposition syndrome
SPECIMEN
WB Whole Blood
RESULT SUMMARY
Negative
RESULT
No reportable variants were detected.
INTERPRETATION
This result decreases the likelihood but does not rule out the presence of a variant in the BAP1 gene. Individuals may have a pathogenic variant in the BAP1 gene that is not detectable by the methods utilized. Additionally, pathogenic variants in other genes not interrogated by this assay may cause a similar phenotype. This result should be interpreted in the context of clinical findings, family history, and other laboratory testing. A genetic consultation may be of benefit.
METHOD
Next generation sequencing (NGS) and/or Sanger sequencing is performed to test for the presence of variants in coding regions and intron/exon boundaries of the genes analyzed, as well as some other regions that have known pathogenic variants. The human genome reference GRCh37/hg19 build was used for sequence read alignment. At least 99% of the bases are covered at a read depth over 30X. Sensitivity is estimated at above 99% for single nucleotide variants, above 94% for indels less than 40 base pairs (bp), above 95% for deletions up to 75 bp and insertions up to 47 bp. NGS, and/or a polymerase chain reaction (PCR)-based quantitative method is performed to test for the presence of deletions and duplications in the genes analyzed. See the Genes Analyzed field for a list of genes tested. There may be regions of genes that cannot be effectively evaluated for sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high GC content, and repetitive sequences. Confirmation of select reportable variants was performed by alternate methodologies based on internal laboratory criteria. See www.mayocliniclabs.com BAP1Z for details regarding genes with regions not routinely covered.
GENES ANALYZED
BAP1
PATIENT NAME TESTRNV, IMPLEMENTATION



PATIENT NAME TESTRNV, IMPLEMENTATION				ORDER NUMBER L002000200
PATIENT ID X100366871	DATE OF BIRTH 02/19/1968	AGE 53 Y	SEX Female	REQUESTED BY CLIENT TESTING
COLLECTED 9/2/2021, 7:00 AM	RECEIVED 9/3/2021, 8:09 AM	REPORTED 9/3/2021, 11:56 AM		
The collected, received, and reported dates and times on the report are in the time zone of the performing location.				CLIENT ORDER NUMBER X100366871
7028846 MCL RochesterCampus Rochester MN 55901				CLIENT MRN 321

DISCLAIMER

Clinical Correlations

If testing was performed because of a clinically significant family history it is often useful to first test an affected family member. Detection of a reportable variant(s) in an affected family member would allow for more informative testing of at risk individuals.

To discuss the availability of further testing options or for assistance in the interpretation of these results, Mayo Clinic Laboratory genetic counselors can be contacted at 1-800-533-1710.

Technical Limitations

Next generation sequencing may not detect all types of genomic variants. In rare cases, false negative or false positive results may occur. The depth of coverage may be variable for some target regions, but assay performance below the minimum acceptable criteria or for failed regions will be noted. Given these limitations, negative results do not rule out the diagnosis of a genetic disorder. If a specific clinical disorder is suspected, evaluation by alternative methods can be considered.

There may be regions of genes that cannot be effectively evaluated for sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high GC content, and repetitive sequences. Confirmation of select reportable variants was performed by alternate methodologies based on internal laboratory criteria.

Additionally, low level mosaic variants may not be detected.

This test is not designed to differentiate between somatic and germline variants. If there is a possibility that any detected variant is somatic, additional testing may be necessary to clarify the significance of results.

Reclassification of Variants Policy

See www.mayocliniclabs.com BAP1Z for information regarding the laboratory's policy for reclassification of variants.

Variant Evaluation

Variant curation is performed using published ACMG-AMP recommendations as a guideline. Other gene-specific guidelines may also be considered. Variants classified as benign or likely benign are not reported.

Results from in silico evaluation tools may change over time and should be interpreted with caution and professional clinical judgment.

TEST CLASSIFICATION

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.



1-800-533-1710
BAP1 Full Gene Analysis

BAP1Z

PATIENT NAME TESTRNV, IMPLEMENTATION				ORDER NUMBER L002000200
PATIENT ID X100366871	DATE OF BIRTH 02/19/1968	AGE 53 Y	SEX Female	REQUESTED BY CLIENT TESTING
COLLECTED 9/2/2021, 7:00 AM	RECEIVED 9/3/2021, 8:09 AM	REPORTED 9/3/2021, 11:56 AM		
The collected, received, and reported dates and times on the report are in the time zone of the performing location. 7028846 MCL RochesterCampus Rochester MN 55901				CLIENT ORDER NUMBER X100366871
				CLIENT MRN 321
RELEASED BY				
Wei Shen, Ph.D.				

Code : MCR Laboratory : Mayo Clinic Laboratories - Rochester Main Campus
Lab Director : WILLIAM G MORICE, II MD, PhD CLIA Certificate : 24D0404292

Address : 200 FIRST STREET SW
ROCHESTER MN 55905