



SAMPLEREPORT, LLTOF ABNORMAL		MEDICAL RECORD # (PATIENT ID) M133691601	ORDER # B004020471
DOB	05/13/1976	CLIENT ID/WARD	7028846
SEX	Male	CLIENT/NAME WARD	DLMP Rochester
CLIENT MRN	0000789	CITY, ST, ZIP	Rochester
REQUESTED BY	CLIENT CLEINT	MN	
		CLIENT ORDER #	0000789
		DATE COLLECTED	11/4/2016 12:23 PM
		DATE RECEIVED	11/4/2016 12:48 PM
		DATE REPORTED	11/4/2016 3:46 PM

Leukemia/Lymphoma; Tech Only Flow

Special Studies:

Peripheral blood, flow cytometric immunophenotyping:

An interpretation was not performed on this case, but will be done by a local pathologist. Care providers seeking assistance on interpretation or would like to request additional testing should call Mayo Medical Laboratories at 1-800-533-1710. Please see associated analysis document.

Quality Assessment: Specimen received within validated guidelines.

A Wright-Giemsa stained slide prepared from the flow cytometry specimen was examined for quality purposes.

Flow cytometry analysis performed with antibodies to the following antigens:

Triage panel: CD3, CD10, CD16, CD19, CD34, CD45 and kappa and lambda surface light chains. B-cell panel: CD5, CD11c, CD19, CD20, CD200, CD22, CD23, CD38, CD45, CD103 and kappa and lambda surface light chains.

Disclaimer:

This test was developed using an analyte specific reagent. Its performance characteristics were determined by Mayo Clinic in a manner consistent with CLIA requirements. This test has not been cleared or approved by the U.S. Food and Drug Administration.

PATIENT NAME: SAMPLEREPORT, LLTOF ABNORMAL

Page 1 of 6

The collected, received, and reported dates and times on the report are in the time zone of the performing location.

MCR Mayo Clinic Laboratories - Rochester Main Campus

200 FIRST STREET SW

ROCHESTER

MN

55905-0001


WILLIAM G MORICE, II, MD, PhD

Patient Name: SAMPLEREPORT, LLTOF...	Mayo Order: 16-YRYH6_LLTOF	Acquisition Date: 04-NOV-2016 11:52:31	Instrument: V33896201643 Ray
Lymphocytes: 70% WBCs	CD3+ T-cells: 15% Lymphs	CD19+ Lymphs 81% Lymphs	Lymphosum: 98
CD34+ Blasts: 0% WBCs	CD3+/CD16+: 0% Lymphs	Lymphs WBCs	Kappa/Lambda
	CD16+ NK-cells: 3% Lymphs	Kappa: 80% 56%	(K/L: 112.38)
		Lambda: 1% 0%	

Triage panel: CD3, CD10, CD16, CD19, CD34, CD45, kappa, and lambda immunoglobulin light chains.


Calculations use Lymphs and Analyzed WBCs for respective denominators.

Kappa/Lambda Ratio will only work if both populations are present.




Friday, November 04, 2016 3:28 PM


Page 1 of 2



Patient Name: SAMPLEREPORT, LL...	Mayo Order: 16-YRYH6_BPAN	Acquisition: 04-NOV-2016 14:12:03	Instrument: V96300273-Barry
Lymphs: 65% of WBCs	CD19+ Lymphs: 86% of Lymphs		

B-cell panel: Tube 1: CD5, CD19, CD20, CD23, CD38, CD45 and kappa and lambda immunoglobulin light chains. Tube 2: CD11c, CD19, CD20, CD22, CD45, CD103 and CD200.





SAMPLEREPORT, LLTOF ABNORMAL

MEDICAL RECORD # (PATIENT ID) M133691601

ORDER # B004020471

DOB 05/13/1976

CLIENT ID/WARD 7028846

CLIENT ORDER # 0000789

SEX Male

CLIENT/NAME WARD DLMP Rochester

DATE COLLECTED 11/4/2016 12:23 PM

CLIENT MRN 0000789


CITY, ST, ZIP Rochester

DATE RECEIVED 11/4/2016 12:48 PM

REQUESTED BY CLIENT CLEINT

MN

DATE REPORTED 11/4/2016 3:46 PM



Page 3 of 3

PATIENT NAME: SAMPLEREPORT, LLTOF ABNORMAL

Page 6 of 6

The collected, received, and reported dates and times on the report are in the time zone of the performing location.