

Reporting Title: JAK2 V617F Mutation Detection, BM
Performing Location: Rochester

Shipping Instructions:
Specimen must arrive within 7 days of collection.

Specimen Requirements:

Container/Tube:
Preferred: Lavender top (EDTA)
Acceptable: Yellow top (ACD)
Specimen Volume: 2 mL
Collection Instructions:

- 1. Invert several times to mix bone marrow.
- 2. Send bone marrow specimen in original tube. **Do not aliquot.**

Forms:

- 1. [Hematopathology Patient Information](#) (T676)
- 2. If not ordering electronically, complete, print, and send a [Hematopathology/Cytogenetics Test Request](#) (T726) with the specimen.

Specimen Type	Temperature	Time	Special Container
Bone Marrow	Ambient (preferred)	7 days	PURPLE OR PINK TOP/EDTA
	Refrigerated	7 days	PURPLE OR PINK TOP/EDTA

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
31157	JAK2 V617F Mutation Detection, BM	Alphanumeric		72333-8
39723	JAK2 Result	Alphanumeric		53761-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81270-JAK2 (Janus kinase 2) (eg, myeloproliferative disorder) gene analysis, p.Val617Phe (V617F) variant

Reference Values:

An interpretive report will be provided.