

# Reporting Title: Chromogranin A, S Performing Location: Rochester

#### **Specimen Requirements:**

Patient Preparation: For at least 2 weeks before specimen collection, patient should stop taking proton pump inhibitor medications. Supplies: Sarstedt Aliquot Tube, 5 mL (T914) **Collection Container/Tube:** Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial Specimen Volume: 0.5 mL

Collection Information: Centrifuge and aliquot serum into plastic vial. Do not submit in original tube.

# Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen: -General Request (T239)

-Oncology Test Request (T729)

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	90 days	
	Ambient	48 hours	
	Refrigerated	48 hours	

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
CGAK	Chromogranin A, S	Numeric	ng/mL	9811-1

LOINC® and CPT codes are provided by the performing laboratory.

#### **Supplemental Report:**

No

# **CPT Code Information:**

86316

# **Reference Values:**

<93 ng/mL Reference values apply to all ages.