

## **Test Definition: PLABN**

Platelet Antibody Screen, Serum

Reporting Title: Platelet Ab Screen, S

Performing Location: Rochester

#### **Ordering Guidance:**

For neonate testing, consider sending a maternal specimen instead of a neonate specimen as unbound platelet antibodies may not be detected in the neonate serum.

This test **is not recommended for** the diagnosis of immune thrombocytopenia or autoimmune thrombocytopenia. Tests that are optimized to detect antibodies bound to the platelets will be useful in these situations; cell-bound platelet antibody (direct) test is strongly recommended.

## **Necessary Information:**

If ordering electronically, answer all prompt questions for timely result reporting:

- 1. Reason for request is required for result interpretation.
- a. Use provided diagnosis options if appropriate. If specific diagnosis is unknown, select the generic answer of alloimmune thrombocytopenia.
- b. Record only the diagnosis pertaining to this test.
- c. Record diagnosis description instead of code.
- 2. Indicate if patient has had intravenous immunoglobulin (IVIg) therapy in the last month: Yes or No

#### **Specimen Requirements:**

Patient Preparation: Do not collect within 72 hours of a platelet transfusion. Transfused platelets will interfere with this

assay.

**Collection Container/Tube:** Red top **Submission Container/Tube:** Plastic vial

Specimen Volume: 1.5 mL

Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

#### Forms:

Platelet Antibody Screen, Serum Patient Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Frozen (preferred)	365 days	
	Refrigerated	48 hours	

## Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
PLABN	PTL15	Reason for request?:	Answer List	Yes
		PTR, Refractory to Platelet Transfusion		
		NAIT, Neonatal Immune Thrombocytopenia		
		PTP, Post-Transfusion Purpura		
		Alloimmune Thrombocytopenia		
		ITP, Order CellBound PLT Ab instead		
		Seconday ITP, order Cell Bound instead		



# **Test Definition: PLABN**

Platelet Antibody Screen, Serum

PLABN	PTL16	IVIg in last month?:	Answer List	Yes
		• Yes		
		• No		

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
PTL01	Overall Result	Alphanumeric		24375-8
PTL02	Interpretation	Alphanumeric		59466-3
PTL03	GPIIb/IIIa (Cell-1)	Alphanumeric		48505-2
PTL04	GPIIb/IIIa (Cell-2)	Alphanumeric		48505-2
PTL05	GPIa/IIa (Cell-1)	Alphanumeric		47084-9
PTL06	GPIa/IIa (Cell-2)	Alphanumeric		47084-9
PTL07	GPIb/IX	Alphanumeric		48506-0
PTL08	GPIV	Alphanumeric		87757-1
PTL09	HLA Class I	Alphanumeric		95269-7
PTL15	Reason for request?	Alphanumeric		29548-5
PTL16	IVIg in last month?	Alphanumeric		95268-9

**LOINC®** and CPT codes are provided by the performing laboratory.

**Supplemental Report:** 

No

**CPT Code Information:** 

86022

**Reference Values:** 

Not applicable