

Test Definition: CLZ

Clozapine, Serum

Reporting Title: Clozapine, S **Performing Location:** Rochester

Ordering Guidance:

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL **Collection Instructions:**

1. Collect specimen immediately before next scheduled dose.

2. Within 2 hours of collection, centrifuge and aliquot serum into a plastic vial.

Additional Information: Therapeutic range (trough level) applies to specimens collected immediately prior to next dose.

Forms:

If not ordering electronically, complete, print, and send a <u>Therapeutics Test Request</u> (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Serum Red	Refrigerated (preferred)	28 days		
	Frozen	28 days		
	Ambient	7 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
33841	Norclozapine	Numeric	ng/mL	10992-6
33842	Clozapine+Norclozapine Total	Alphanumeric	ng/mL	12375-2
33840	Clozapine	Numeric	ng/mL	6896-5

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80159

Reference Values:

Clozapine

Therapeutic range: 350-600 ng/mL

Norclozapine



Test Definition: CLZ

Clozapine, Serum

Therapeutic range: Not well established

Clozapine + Norclozapine

Therapeutic range: Not well established