

Reporting Title: ANNA-1 Titer, S
Performing Location: Rochester

Specimen Requirements:

- Only orderable as a reflex. For more information see:
- AIAES / Axonal Neuropathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - DMS2 / Dementia, Autoimmune/Paraneoplastic Evaluation, Serum
 - DYS2 / Dysautonomia, Autoimmune/Paraneoplastic Evaluation, Serum
 - ENS2 / Encephalopathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - EPS2 / Epilepsy, Autoimmune/Paraneoplastic Evaluation, Serum
 - GID2 / Gastrointestinal Dysmotility, Autoimmune/Paraneoplastic Evaluation, Serum
 - MAS1 / Myelopathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - MDS2 / Movement Disorder, Autoimmune/Paraneoplastic Evaluation, Serum
 - PAVAL / Paraneoplastic, Autoantibody Evaluation, Serum
 - PCDES / Pediatric Autoimmune Encephalopathy/CNS Disorder Evaluation, Serum

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
43431	ANNA-1 Titer, S	Alphanumeric	titer	94342-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
86256

Reference Values:

- Only orderable as a reflex. For more information see:
- AIAES / Axonal Neuropathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - DMS2 / Dementia, Autoimmune/Paraneoplastic Evaluation, Serum
 - DYS2 / Dysautonomia, Autoimmune/Paraneoplastic Evaluation, Serum
 - ENS2 / Encephalopathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - EPS2 / Epilepsy, Autoimmune/Paraneoplastic Evaluation, Serum
 - GID2 / Gastrointestinal Dysmotility, Autoimmune/Paraneoplastic Evaluation, Serum

-
- MAS1 / Myelopathy, Autoimmune/Paraneoplastic Evaluation, Serum
 - MDS2 / Movement Disorder, Autoimmune/Paraneoplastic Evaluation, Serum
 - PAVAL / Paraneoplastic, Autoantibody Evaluation, Serum
 - PCDES / Pediatric Autoimmune Encephalopathy/CNS Disorder Evaluation, Serum

<1:240

Neuron-restricted patterns of IgG staining that do not fulfill criteria for antineuronal nuclear antibody type 1 may be reported as "unclassified anti-neuronal IgG." Complex patterns that include nonneuronal elements may be reported as "uninterpretable."