

Circulating Immune Complex, C1q Binding

Reporting Title: CIC, C1q Binding

Performing Location: ARUP Laboratories

Specimen Requirements:

Draw blood in a plain, red-top tube(s). (Serum gel tube is acceptable.) Separate from cells and freeze immediately. Send 1 mL of serum frozen in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z2636	CIC, C1q Binding	Alphanumeric		27831-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86332

Reference Values:

0.0-3.9 ug Eq/mL

Less than or equal to 3.9 ug Eq/mL is considered negative for circulating complement binding immune complexes.

Circulating immune complexes may be found without any evident pathology and positive results do not necessarily implicate the immune complex in a disease process.