

Reporting Title: Pantothenic Acid (B-5) Bioassay

Performing Location: BioAgilytix Diagnostics

Specimen Requirements:

Specimen Type: Serum

Container/Tube: SST

Specimen Volume: 2 mL

Collection Instructions: Draw blood in SST (serum separator tube). Allow sample to clot for 30 minutes at room temperature. Spin down and send 1 mL serum frozen in an amber vial (T915) to protect from light.

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	21 days	LIGHT PROTECTED
	Refrigerated	7 days	LIGHT PROTECTED

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FPAB	Pantothenic Acid (B-5) Bioassay	Alphanumeric	ug/L	2722-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84591

Reference Values:

Adult Reference Range

>10 Years 37-147 ug/L

Pediatric Reference Range

< or =1 Year 3.45 to 825 ug/L

>1 year to 10 Years 3.45 to 229.2 ug/L