

Reporting Title: Annatto Seed IgE
Performing Location: Eurofins Viracor

Specimen Requirements:

Container/Tube: Red or SST

Preferred: Red top tube

Acceptable: Serum gel tube

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3074	Annatto Seed IgE	Alphanumeric		50638-6
Z3075	Class	Alphanumeric		103035-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L