

Reporting Title: IL-8, Serum
Performing Location: Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum
Collection Container/Tube: Red or SST
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FIL8S	IL-8, Serum	Alphanumeric		33211-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83520

Reference Values:

<12.4 pg/mL