

**Reporting Title:** IL-4, Serum  
**Performing Location:** Eurofins Viracor

**Specimen Requirements:**

**Specimen Type:** Serum  
**Collection Container/Tube:** Red or SST  
**Submission Container/Tube:** Plastic vial  
**Specimen Volume:** 1 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum frozen in a plastic vial.

| Specimen Type | Temperature | Time    | Special Container |
|---------------|-------------|---------|-------------------|
| Serum         | Frozen      | 14 days |                   |

**Result Codes:**

| Result ID | Reporting Name | Type         | Unit | LOINC®  |
|-----------|----------------|--------------|------|---------|
| FIL4S     | IL-4, Serum    | Alphanumeric |      | 27161-9 |

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83520

**Reference Values:**

<0.5 pg/mL