

Banana IgG

# **Reporting Title:** Banana IgG **Performing Location:** Eurofins Viracor

## **Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red Top

Specimen Volume: 0.5 mL

**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FBANG	Banana IgG	Alphanumeric		60352-2

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

## **CPT Code Information:**

86001

#### **Reference Values:**

Reference ranges have not been established for food-specific IgG tests.