

# **Reporting Title:** Blueberry IgG **Performing Location:** Eurofins Viracor

### **Specimen Requirements:**

Specimen Type: Serum

Container/Tube: Red or SST Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain, red-top tube(s) or serum gel tube(s). Spin down and send 0.5 mL of serum refrigerate in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
FBLUG	Blueberry IgG	Alphanumeric		58035-7

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

# **CPT Code Information:**

86001

## **Reference Values:**

Reference ranges have not been established for food-specific IgG tests.