

Test Definition: FAVCG

Avocado IgG

Reporting Title: Avocado IgG **Performing Location:** Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum Container/Tube: Red Top Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5

mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FAVCG	Avocado IgG	Alphanumeric		60398-5

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86001

Reference Values:

Reference ranges have not been established for food-specific IgG tests.