

Reporting Title: Flounder IgE
Performing Location: Eurofins Viracor

Specimen Requirements:
Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3815	Flounder IgE	Alphanumeric		7316-3
Z3816	CLASS	Alphanumeric		102311-8

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
86003

Reference Values:
<0.35 kU/L