

Reporting Title: Gelatin Porcine IgE
Performing Location: Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3876	Gelatin Porcine IgE	Alphanumeric		16103-4
Z3877	CLASS	Alphanumeric		102745-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L