

Reporting Title: Bird Fancier's Precipitin Panel I

Performing Location: Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum

Container/Tube: Red or SST

Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3892	Parakeet Droppings Gel Diffusion	Alphanumeric		24503-5
Z3894	Parrott Droppings Gel Diffusion	Alphanumeric		24501-9
Z3890	Cockatiel Droppings Gel Diffusion	Alphanumeric		23843-6
Z3891	Finch Droppings Gel Diffusion	Alphanumeric		Not Provided
Z3893	Parakeet Serum Gel Diffusion	Alphanumeric		24502-7
Z3888	Canary Droppings Gel Diffusion	Alphanumeric		Not Provided
Z3889	Chicken Serum Gel Diffusion	Alphanumeric		24498-8
Z3895	Parrott Serum Gel Diffusion	Alphanumeric		24500-1
Z3897	Pigeon/Dove Serum Gel Diffusion	Alphanumeric		6733-0
Z3896	Pigeon/Dove Droppings Gel Diffusion	Alphanumeric		24505-0

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86331 x 10

Reference Values:

Negative