

Reporting Title: Aspen (Populus tremuloides) IgE

Performing Location: Eurofins Viracor

Specimen Requirements:

Specimen Type: Serum

Container/Tube: Red Top

Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |
|               | Ambient                  | 28 days  |                   |

Result Codes:

| Result ID | Reporting Name                  | Type         | Unit | LOINC®   |
|-----------|---------------------------------|--------------|------|----------|
| Z4012     | Aspen (Populus tremuloides) IgE | Alphanumeric |      | 7101-9   |
| Z4013     | CLASS                           | Alphanumeric |      | 102152-6 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L