

**Reporting Title:** Donath Landsteiner  
**Performing Location:** ARUP Laboratories

**Specimen Requirements:**

**Specimen Type:** Serum  
**Container/Tube:** Red  
**Specimen Volume:** 3 mL

Draw blood in a plain red-top tube. Maintain specimen at 37 degrees C until serum is separated from cells. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	7 days	
	Ambient	72 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
FDLS	Donath Landsteiner	Alphanumeric		18287-3

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86940, 86941

**Reference Values:**

Negative