

# **Test Definition: PF199**

Fluid

Carbohydrate Antigen 19-9 (CA 19-9), Pleural

Reporting Title: CA 19-9, Pleural Fluid

Performing Location: Rochester

## **Specimen Requirements:**

**Patient Preparation: For 12 hours before specimen collection do not** take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.

Container/Tube: Plain, plastic, screw-top tube

Specimen Volume: 2 mL

#### Forms:

If not ordering electronically, complete, print, and send an Oncology Test Request (T729) with the specimen.

| Specimen Type | Temperature        | Time    | Special Container |  |
|---------------|--------------------|---------|-------------------|--|
| Pleural Fluid | Frozen (preferred) | 90 days |                   |  |
|               | Refrigerated       | 14 days |                   |  |
|               | Ambient            | 7 days  |                   |  |

#### **Result Codes:**

| Result ID | Reporting Name         | Туре    | Unit | LOINC®  |
|-----------|------------------------|---------|------|---------|
| P199      | CA 19-9, Pleural Fluid | Numeric | U/mL | 19163-5 |
| SITE8     | Site                   | Numeric |      | 39111-0 |

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

86301

### **Reference Values:**

An interpretive report will be provided.