

# **Test Definition: SHUR**

Schistosoma Exam, Random, Urine

Reporting Title: Schistosoma Exam, U

Performing Location: Rochester

**Specimen Requirements:** 

Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Clean, plastic urine collection container

Submission Container/Tube: Plastic, 10-mL urine tube

**Specimen Volume:** 10 mL **Collection Instructions:** 

1. Collect a random urine specimen. Preferred time of collection between the hours of 12 noon and 3 p.m. but not

required. A 24-hour urine collection is also acceptable.

2. No preservative.

| Specimen Type | Temperature  | Time   | Special Container |
|---------------|--------------|--------|-------------------|
| Urine         | Refrigerated | 7 days |                   |

## Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description     | Туре       | Reportable |
|---------|-------------|-----------------|------------|------------|
| SHUR    | Q00M0052    | Specimen Source | Plain Text | Yes        |

## **Result Codes:**

| Result ID | Reporting Name      | Туре         | Unit | LOINC®  |
|-----------|---------------------|--------------|------|---------|
| SHUR      | Schistosoma Exam, U | Alphanumeric |      | 10715-1 |

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

87210 87015

## **Reference Values:**

Negative

If positive, organism identified