

Test Definition: VITAP

Vitamin A, Serum

Reporting Title: Vitamin A, S **Performing Location:** Rochester

Specimen Requirements:

Only orderable as part of a profile. For more information see VITAE / Vitamin A and Vitamin E, Serum.

Patient Preparation: Fasting overnight (12-14 hours) (infants: collect specimen prior to next feeding)

Supplies: Amber Frosted Tube, 5 mL (T915)

Collection Container/Tube:

Preferred: Red top **Acceptable:** Serum gel

Submission Container/Tube: Amber vial

Specimen Volume: 0.5 mL

Collection Instructions: Within 24 hours of collection, centrifuge, and aliquot serum into an amber vial to protect from

light.

Specimen Type	Temperature	Time	Special Container	
Serum	Refrigerated (preferred)	28 days	LIGHT PROTECTED	
	Frozen	28 days	LIGHT PROTECTED	
	Ambient	7 days	LIGHT PROTECTED	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
605124	Vitamin A	Numeric	mcg/dL	2923-1

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84590

Reference Values:

Only orderable as part of a profile. For more information see VITAE / Vitamin A and Vitamin E, Serum.

0-6 years: 11.3-64.7 mcg/dL 7-12 years: 12.8-81.2 mcg/dL 13-17 years: 14.4-97.7 mcg/dL > or =18 years: 32.5-78.0 mcg/dL