

**Reporting Title:** N-Methylhistamine, 24 Hr, U

**Performing Location:** Rochester

**Ordering Guidance:**

A random urine collection is preferred for patients with episodic symptoms, for example in the context of allergic reactions, brought on by specific environmental factors. See NMHR / N-Methylhistamine, Random, Urine.

**Specimen Requirements:**

Only orderable as part of a profile. For more information see NMH24 / N-Methylhistamine, 24 Hour, Urine

**Supplies:** Sarstedt 5 mL Aliquot Tube (T914)

**Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 5 mL

**Collection Instructions:**

1. Collect urine for 24 hours.
2. No preservative.

**Urine Preservative Collection Options:**

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
605159	N-Methylhistamine, 24 Hr	Alphanumeric	mcg/g Cr	44340-8

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

82542

**Reference Values:**

Only orderable as part of a profile. For more information see:

NMH24 / N-Methylhistamine, 24 Hour, Urine

0-5 years: 120-510 mcg/g creatinine

6-16 years: 70-330 mcg/g creatinine

---

>16 years: 30-200 mcg/g creatinine