

Reporting Title: TP53 Pre-Analysis Cell Sorting, V
Performing Location: Rochester

Specimen Requirements:
Only orderable as a reflex. For more information see P53CA / Hematologic Neoplasms, *TP53* Somatic Mutation, DNA Sequencing Exons 4-9, Varies.

Specimen Type: Blood
Container/Tube: Lavender top (EDTA) or yellow top (ACD solution B)
Specimen Volume: 3 mL
Collection Instructions:
1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**
3. Label specimen as blood.
Specimen Stability Information: Ambient/Refrigerate <10 days

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	10 days	
	Refrigerated	10 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
607678	TP53 Pre-Analysis Cell Sort	Alphanumeric		No LOINC Needed
607687	Final Diagnosis	Alphanumeric		22637-3

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
88184-Flow cytometry, first cell surface, cytoplasmic or nuclear marker
88185 x 4-Each additional marker

Reference Values:
Only orderable as a reflex. For more information see P53CA / Hematologic Neoplasms, *TP53* Somatic Mutation, DNA Sequencing Exons 4-9, Varies.