

# **Test Definition: CSMRT**

mSMART Plasma Cell Proliferative Disorder, Pre-Analysis Cell Sorting, Bone Marrow

Reporting Title: MPCDS Pre-Analysis Cell Sorting, BM

Performing Location: Rochester

## **Specimen Requirements:**

Only orderable as a reflex. For more information see MSMRT / Mayo Algorithmic Approach for Stratification of Myeloma and Risk-Adapted Therapy Report, Bone Marrow or MSMRD / Myeloma Stratification and Risk-Adapted Therapy with Reflex to Minimal Residual Disease, Bone Marrow

**Specimen Type:** Redirected bone marrow **Preferred:** Yellow top (ACD solution A or B)

Acceptable: Lavender top (EDTA) or green top (heparin)

Specimen Volume: 4 mL

Specimen Type	Temperature	Time	Special Container
Bone Marrow	Ambient (preferred)	4 days	
	Refrigerated	4 days	

# **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
607682	MPCDS Pre-Analysis Cell Sort	Alphanumeric		No LOINC Needed
607688	Final Diagnosis	Alphanumeric		No LOINC Needed

LOINC® and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

#### **CPT Code Information:**

88184-Flow Cytometry; first cell surface, cytoplasmic or nuclear marker 88185 x 5-Flow Cytometry, additional cell surface, cytoplasmic or nuclear marker (each)

### **Reference Values:**

Only orderable as a reflex. For more information see MSMRT / Mayo Algorithmic Approach for Stratification of Myeloma and Risk-Adapted Therapy Report, Bone Marrow or MSMRD / Myeloma Stratification and Risk-Adapted Therapy with Reflex to Minimal Residual Disease, Bone Marrow

An interpretive report will be provided.