

Test Definition: ALDEF

Alpha Defensin, Lateral Flow Assay, Synovial

Fluid

Reporting Title: Alpha Defensin, Synovial Fluid, LFA

Performing Location: Rochester

Ordering Guidance:

This test should **only be used** for patients with a total joint prosthesis.

Necessary Information:

Specimen source and type of joint are required. If not obtained from a prosthetic joint, testing will be canceled.

If ordering electronically, answers must be provided for the order entry questions.

If not ordering electronically, specimen source and type of joint must be provided on the request form.

Specimen Requirements:

Collection Container/Tube: Plain red-top tube **Submission Container/Tube:** Plastic vial

Specimen Volume: 0.6 mL

Specimen Type	Temperature	Time	Special Container
Synovial Fluid	Refrigerated	7 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
ALDEF	JOINT	Is this for a prosthetic joint?:	Answer List	Yes
		• Yes		
		• No		
		Unknown		
ALDEF	ADSRC	Alpha Defensin Source:	Answer List	Yes
		Synovial Fluid, Left Ankle		
		Synovial Fluid, Right Ankle		
		Synovial Fluid, Left Elbow		
		Synovial Fluid, Right Elbow		
		Synovial Fluid, Left Hip		
		Synovial Fluid, Right Hip		
		Synovial Fluid, Left Knee		
		Synovial Fluid, Right Knee		
		Synovial Fluid, Left Shoulder		
		Synovial Fluid, Right Shoulder		
		Synovial Fluid, Left Wrist		
		Synovial Fluid, Right Wrist		
		Unknown		



Test Definition: ALDEF

Alpha Defensin, Lateral Flow Assay, Synovial

Fluid

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
607711	Alpha Defensin, Synovial Fl, LFA	Alphanumeric		94718-4
JOINT	Is this for a prosthetic joint?	Alphanumeric		86955-2
ADSRC	Alpha Defensin Source	Alphanumeric		72923-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83518

Reference Values:

Negative

Reference values apply to all ages.