
Reporting Title: CYP2C19 Genotype, V**Performing Location:** Rochester**Ordering Guidance:**

Testing is available as the single gene assay (this test) and as a part of a psychotropic or focused pharmacogenomics panel.

If multiple pharmacogenomic genotype testing is needed, consider PGXQP / Focused Pharmacogenomics Panel, Varies.

If genotype testing for psychotropic medications is requested, order PSYQP / Psychotropic Pharmacogenomics Gene Panel, Varies.

Specimen Requirements:

Multiple genotype tests can be performed on a single specimen after a single extraction. See [Multiple Genotype Test List](#) for a list of tests that can be ordered together.

Submit only 1 of the following specimens:**Specimen Type:** Whole blood**Container/Tube:** Lavender top (EDTA)**Specimen Volume:** 3 mL**Collection Instructions:**

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not aliquot.**

Specimen Stability Information: Ambient (preferred) 9 days/Refrigerated 30 days**Specimen Type:** Saliva**Patient Preparation:** Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.**Supplies:** Saliva Swab Collection Kit (T786)**Specimen Volume:** 1 Swab**Collection Instructions:** Collect and send specimen per kit instructions.**Specimen Stability Information:** Ambient 30 days**Specimen Type:** Extracted DNA**Container/Tube:** 2 mL screw top tube**Specimen Volume:** 100 mcL (microliters)**Collection Instructions:**

1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.
2. Provide concentration of DNA and volume on tube.

Specimen Stability Information: Frozen (preferred)/Ambient/Refrigerated**Forms:**

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:

-[Informed Consent for Genetic Testing](#) (T576)

-[Informed Consent for Genetic Testing-Spanish](#) (T826)

2. If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

-[Neurology Specialty Testing Client Test Request](#) (T732)

-[Therapeutics Test Request](#) (T831)

-[Cardiovascular Test Request](#) (T724)

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
610089	CYP2C19 Genotype	Alphanumeric		57132-3
610090	CYP2C19 Phenotype	Alphanumeric		79714-2
610567	CYP2C19 Activity Score	Alphanumeric		In Process
610091	Interpretation	Alphanumeric		69047-9
610092	Additional Information	Alphanumeric		48767-8
610093	Method	Alphanumeric		85069-3
610094	Disclaimer	Alphanumeric		62364-5
610095	Reviewed by	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81225

Reference Values:

An interpretive report will be provided.