

**Reporting Title:** HLA-B 5701 Genotype, V  
**Performing Location:** Rochester

**Specimen Requirements:**  
[Multiple genotype tests can be performed on a single specimen after a single extraction. See Multiple Genotype Test List](#)  
for a list of tests that can be ordered together.

**Submit only 1 of the following specimens:**

**Specimen Type:** Whole blood  
**Container/Tube:** Lavender top (EDTA)  
**Specimen Volume:** 3 mL  
**Collection Instructions:**  
1. Invert several times to mix blood.  
2. Send whole blood specimen in original tube. **Do not aliquot.**  
**Specimen Stability Information:** Ambient (preferred) 9 days/Refrigerated 30 days

**Specimen Type:** Saliva  
**Supplies:** Saliva Swab Collection Kit (T786)  
**Patient Preparation:** Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.  
**Specimen Volume:** 1 swab  
**Collection Instructions:** Collect and send specimen per kit instructions.  
**Specimen Stability Information:** Ambient 30 days

**Specimen Type:** Extracted DNA  
**Container/Tube:** 2 mL screw top tube  
**Specimen Volume:** 100 mcL (microliters)  
**Collection Instructions:**  
1. The preferred volume is 100 mcL at a concentration of 50 ng/mcL.  
2. Include concentration and volume on tube.  
**Specimen Stability Information:** Frozen (preferred) 1 year/Ambient/Refrigerated

**Forms:**  
1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:  
[-Informed Consent for Genetic Testing \(T576\)](#)  
[-Informed Consent for Genetic Testing-Spanish \(T826\)](#)  
2. [If not ordering electronically, complete, print, and send a Therapeutics Test Request \(T831\)](#) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
610672	HLA-B *57:01 Genotype	Alphanumeric		50956-2
610673	HLA-B *57:01 Phenotype	Alphanumeric		93308-5
610674	Interpretation	Alphanumeric		69047-9
610675	Additional Information	Alphanumeric		48767-8
610676	Method	Alphanumeric		85069-3
610677	Disclaimer	Alphanumeric		62364-5
610678	Reviewed by	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81381

Reference Values:

Negative

An interpretive report will be provided.