

**Reporting Title:** Creatinine, Random, U  
**Performing Location:** Rochester

- Specimen Requirements:**  
**Supplies:** Urine Tubes, 10 mL (T068)  
**Container/Tube:** Plastic, 10-mL urine tube  
**Specimen Volume:** 6 mL  
**Collection Instructions:**  
1. Collect a random urine specimen.  
2. No preservative.

**Forms:**  
If not ordering electronically, complete, print, and send a [Renal Diagnostics Test Request](#) (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	30 days	
	Frozen	30 days	
	Ambient	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
RCTUR	Creatinine, Random, U	Numeric	mg/dL	2161-8

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
82570

**Reference Values:**  
> or =18 years old: 16-326 mg/dL

Reference values have not been established for patients who are <18 years of age.