

Reporting Title: Chronic Hepatitis B Screen, S

Performing Location: Rochester

Necessary Information:

Date of collection is required.

Specimen Requirements:

Patient Preparation: For 24 hours before specimen collection, patient **should not** take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 1.2 mL

Collection Instructions:

- Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
- Aliquot serum into plastic vial.

Forms:

If not ordering electronically, complete, print, and send 1 of the following:

[-Gastroenterology and Hepatology Test Request](#) (T728)

[-Infectious Disease Serology Test Request](#) (T916)

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	90 days	
	Refrigerated	6 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
HBAGS	HBs Antigen Scrn, S	Alphanumeric		5196-1

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Components:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
HBGSN	HBs Antigen Scrn, S	1	87340	Yes	Yes

CPT Code Information:

87340
G0499-(if appropriate)
87350 (if appropriate)
87341 (if appropriate)
86707 (if appropriate)

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
EAG	Hepatitis Be Ag, S	1	87350	No	Yes
HBGSC	HBs Antigen Screen Confirmation, S	1	87341	No	No
HEAB	HBe Antibody, S	1	86707	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
EAG	EAG	Hepatitis Be Ag, S	Alphanumeric		13954-3
HEAB	HEAB	HBe Antibody, S	Alphanumeric		33463-1
HBGSC	HBGSC	HBs Antigen Screen Confirmation, S	Alphanumeric		7905-3

Reference Values:

Negative
See [Viral Hepatitis Serologic Profiles](#)