

Reporting Title: CA 19-9, Peritoneal Fluid
Performing Location: Rochester

Specimen Requirements:
Patient Preparation: For 12 hours before specimen collection do not take multivitamins or dietary supplements containing biotin (vitamin B7), which is commonly found in hair, skin, and nail supplements and multivitamins.
Container/Tube: Plain, plastic, screw-top tube
Specimen Volume: 2 mL

Forms:
[If not ordering electronically, complete, print, and send an Oncology Test Request](#) (T729) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------|---------|-------------------|
| Peritoneal | Frozen (preferred) | 90 days | |
| | Ambient | 7 days | |
| | Refrigerated | 7 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|---------------------------|---------|------|---------|
| 199PN | CA 19-9, Peritoneal Fluid | Numeric | U/mL | 50781-4 |
| SITEE | Site | Numeric | | 39111-0 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
86301

Reference Values:
An interpretive report will be provided.