

## **Test Definition: DWPAN**

Comprehensive Distal Weakness Gene Panel, Varies

# **Reporting Title:** Distal Weakness Gene Panel **Performing Location:** Rochester

#### Ordering Guidance:

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

Customization of this panel and single gene analysis for any gene present on this panel are available. For more information see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

#### Shipping Instructions:

Specimen preferred to arrive within 96 hours of collection.

#### **Specimen Requirements:**

**Patient Preparation:** A previous bone marrow transplant from an allogenic donor will interfere with testing. Call 800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Specimen Type: Whole blood

Container/Tube: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant

Specimen Volume: 3 mL

**Collection Instructions:** 

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred)/Refrigerated

#### Forms:

1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file.

The following documents are available:

-Informed Consent for Genetic Testing (T576)

-Informed Consent for Genetic Testing (Spanish) (T826)

2. Molecular Genetics: Neurology Patient Information

3. If not ordering electronically, complete, print, and send a <u>Neurology Specialty Testing Client Test Request</u> (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

#### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
617546	Test Description	Alphanumeric		62364-5
617547	Specimen	Alphanumeric		31208-2



## **Test Definition: DWPAN**

### Comprehensive Distal Weakness Gene Panel, Varies

617548	Source	Alphanumeric	31208-2
617549	Result Summary	Alphanumeric	50397-9
617550	Result	Alphanumeric	82939-0
617551	Interpretation	Alphanumeric	69047-9
618178	Additional Results	Alphanumeric	82939-0
617552	Resources	Alphanumeric	99622-3
617553	Additional Information	Alphanumeric	48767-8
617554	Method	Alphanumeric	85069-3
617555	Genes Analyzed	Alphanumeric	48018-6
617556	Disclaimer	Alphanumeric	62364-5
617557	Released By	Alphanumeric	18771-6

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

Supplemental

## **CPT Code Information:** 81443

81443

#### **Reference Values:**

An interpretive report will be provided.