

Reporting Title: Hazelnut-Food Component Reflex, S
Performing Location: Rochester

Ordering Guidance:
For a listing of allergens available for testing, see [Allergens - Immunoglobulin E \(IgE\) Antibodies](#).

Specimen Requirements:
Collection Container/Tube:
Preferred: Serum gel
Acceptable: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

Forms:
If not ordering electronically, complete, print, and send an [Allergen Test Request](#) (T236) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	90 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
NUTH1	Hazelnut-Food, IgE, S	Numeric	kU/L	6136-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
86003

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
NUTHX	Hazelnut-Food Components, IgE, S	1	86003	No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
---------	-----------	----------------	------	------	--------

Test Definition: NUTHR

Hazelnut-Food, IgE with Reflex to
Hazelnut-Food Components, IgE, Serum

NUTHX	INUTH	Hazelnut-Food IgE Ab Interpretation	Alphanumeric		69048-7
NUTHX	A1COR	Cor a 1 (Hazelnut-Food), IgE, S	Numeric	kU/L	69421-6
NUTHX	A8COR	Cor a 8 (Hazelnut-Food), IgE, S	Numeric	kU/L	58753-5
NUTHX	A9COR	Cor a 9 (Hazelnut-Food), IgE, S	Numeric	kU/L	65765-0
NUTHX	A14CO	Cor a 14 (Hazelnut-Food), IgE, S	Numeric	kU/L	81788-2

Reference Values:

Class	IgE kU/L	Interpretation
0	<0.10	Negative
0/1	0.10-0.34	Borderline/Equivocal
1	0.35-0.69	Equivocal
2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive
6	> or =100	Strongly positive

Concentrations of 0.70 kU/L or more (class 2 and above) will flag as abnormally high.
Reference values apply to all ages.