

Reporting Title: Hazelnut-Food Components, IgE, S
Performing Location: Rochester

Specimen Requirements:
Only orderable as a reflex. For more information see NUTHR / Hazelnut-Food, IgE with Reflex to Hazelnut-Food Components, IgE, Serum.

Container/Tube:
Preferred: Serum gel
Acceptable: Red top
Submission Container/Tube: Plastic vial
Specimen Volume: 1 mL
Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	90 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
INUTH	Hazelnut-Food IgE Ab Interpretation	Alphanumeric		69048-7
A1COR	Cor a 1 (Hazelnut-Food), IgE, S	Numeric	kU/L	69421-6
A8COR	Cor a 8 (Hazelnut-Food), IgE, S	Numeric	kU/L	58753-5
A9COR	Cor a 9 (Hazelnut-Food), IgE, S	Numeric	kU/L	65765-0
A14CO	Cor a 14 (Hazelnut-Food), IgE, S	Numeric	kU/L	81788-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
86003
86008 x 3

Reference Values:
Only orderable as a reflex. For more information see NUTHR / Hazelnut-Food, IgE with Reflex to Hazelnut-Food Components, IgE, Serum.

Class	IgE kU/L	Interpretation
0	<0.10	Negative
0/1	0.10-0.34	Borderline / Equivocal

1	0.35-0.69	Equivocal
2	0.70-3.49	Positive
3	3.50-17.4	Positive
4	17.5-49.9	Strongly positive
5	50.0-99.9	Strongly positive
6	> or =100	Strongly positive

Concentrations of 0.70 kU/L or more (class 2 and above) will flag as abnormally high.
Reference values apply to all ages.