

# **Test Definition: HGEMS**

Hydroxyglutaric Acids, Glutaric Acid, Ethylmalonic Acid, and Methylsuccinic Acid, Serum

#### **Reporting Title:** HGEM, S **Performing Location:** Rochester

## **Necessary Information:**

Patient's age is required.

### **Specimen Requirements:**

Collection Container/Tube: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.1 mL Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

#### Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	90 days	
	Frozen	90 days	
	Ambient	4 days	

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
36055	2-OH Glutaric acid	Numeric	nmol/mL	69845-6
36056	3-OH Glutaric acid	Numeric	nmol/mL	69851-4
36057	Glutaric acid	Numeric	nmol/mL	27301-1
36058	Methylsuccinic acid	Numeric	nmol/mL	69829-0
36059	Ethylmalonic acid	Numeric	nmol/mL	79476-8
36060	Interpretation (HGEMS)	Alphanumeric		59462-2
36061	Reviewed By	Alphanumeric		18771-6

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

## Supplemental Report:

No

# **CPT Code Information:** 83918

#### **Reference Values:**

2-OH Glutaric acid: < or =4.5 nmol/mL



# **Test Definition: HGEMS**

Hydroxyglutaric Acids, Glutaric Acid, Ethylmalonic Acid, and Methylsuccinic Acid, Serum

3-OH Glutaric acid: < or =0.7 nmol/mL Glutaric acid: < or =0.8 nmol/mL Methylsuccinic acid: < or =0.3 nmol/mL Ethylmalonic acid: < or =1.5 nmol/mL