

Reporting Title: Oxycodone Screen, U

Performing Location: Rochester

**Ordering Guidance:**

For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order OXYSX / Oxycodone Screen, Chain of Custody, Random, Urine.

Additional drug panels and specific requests are available; call 800-533-1710 or 507-266-5700.

**Additional Testing Requirements:**

If urine creatinine is required or adulteration of the sample is suspected, the following test should also be ordered, ADULT / Adulterants Survey, Random, Urine.

**Specimen Requirements:**

**Supplies:** Sarstedt 5 mL Aliquot Tube (T914)

**Collection Container/Tube:** Plastic urine container

**Submission Container/Tube:** Plastic, 5-mL tube

**Specimen Volume:** 2 mL

**Collection Instructions:**

1. Collect a random urine specimen.
2. No preservative.

**Forms:**

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
62623	Oxycodone	Alphanumeric	ng/mL	19642-8

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

80307

**Reflex Tests:**

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
OXYCU	Oxycodone w/metabolite Conf, U	1	80365	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
OXYCU	62616	Oxycodone-by LC-MS/MS	Alphanumeric	ng/mL	16249-5
OXYCU	35965	Oxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	17395-5
OXYCU	36022	Oxycodone Interpretation	Alphanumeric		69050-3
OXYCU	42010	Noroxycodone-by LC-MS/MS	Alphanumeric	ng/mL	61425-5
OXYCU	42011	Noroxymorphone-by LC-MS/MS	Alphanumeric	ng/mL	90894-7

Reference Values:

Negative

Screening cutoff concentration:

Oxycodone: 100 ng/mL