

**Reporting Title:** Chromogenic FIX, P  
**Performing Location:** Rochester

**Ordering Guidance:**  
Coagulation testing is highly complex, often requiring the performance of multiple assays and correlation with clinical information. For that reason, ALBLD / Bleeding Diathesis Profile, Limited, Plasma is recommended.

**Necessary Information:**

**Specimen Requirements:**  
**Specimen Type:** Platelet-poor plasma  
**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)  
**Submission Container/Tube:** Polypropylene vial preferred  
**Specimen Volume:** 1 mL

**Collection Instructions:**

- For complete instructions, see [Coagulation Guidelines for Specimen Handling and Processing](#)
- Centrifuge, transfer all plasma into a plastic vial, and centrifuge plasma again.
- Aliquot plasma into a separate plastic vial leaving 0.25 mL in the bottom of centrifuged vial.
- Freeze plasma immediately (no longer than 4 hours after collection) at -20 degrees C or, ideally, at -40 degrees C or below.

**Additional Information:**

- Double-centrifuged specimen is critical for accurate results as platelet contamination may cause spurious results.
- Each coagulation assay requested should have its own vial.

**Forms:**  
[If not ordering electronically, complete, print, and send a Coagulation Test Request](#) (T753) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
CH9	Chromogenic FIX, P	Numeric	%	88449-4

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
85130

---

**Reference Values:**

65-140%

Chromogenic factor IX activity generally correlates with the one-stage FIX activity. In full term/premature neonates, infants, children, and adolescents the one-stage FIX activity\* is similar to adults. However, no similar data for chromogenic FIX activity are available.(Appel IM, Grimminck B, Geerts J, Stigter R, Cnossen MH, Beishuizen A. Age dependency of coagulation parameters during childhood and puberty. J Thromb Haemost. 2012;10(11):2254-2263)

\*See Pediatric Hemostasis References section in [Coagulation Guidelines for Specimen Handling and Processing](#)