

# **Test Definition: LEPDT**

Leptospira, IgM, Serum

Reporting Title: Leptospira, IgM, S

Performing Location: Rochester

**Specimen Requirements:** Collection Container/Tube:

**Preferred:** Serum gel **Acceptable:** Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 0.3 mL

**Collection Instructions:** Centrifuge and aliquot serum into plastic vial.

Additional Information: If acute and convalescent specimens are being obtained to determine seroconversion, they

should be collected 2 or more weeks apart.

#### Forms:

If not ordering electronically, complete, print, and send <u>Infectious Disease Serology Test Request</u> (T916) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	30 days	
	Frozen	30 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
65183	Leptospira, IgM, S	Alphanumeric		23201-7

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

## **CPT Code Information:**

86720

## **Reference Values:**

Negative